


Ampullenspritze Cartridge Syringes

Available in European and American Both Threads



MSI-1743-10
Flat Wing Non-Folding Type
Anesthesia syringe. Cylinder
without Aspiration Metric / Us
Thread - 1.8 CC & 2.2 CC



MSI-1743-11
Flat Wing Folding Type
Anesthesia syringe. Cylinder
without Aspiration Metric / Us
Thread - 1.8 CC & 2.2 CC



MSI-1743-12
Flat Wing Cartridge Syringe
for Anesthesia, folding type. Cylinder
without Aspiration Metric / Us
Thread - 1.8 CC & 2.2 CC



MSI-1743-13
Flat Wing Folding Type
Anesthesia syringe. Cylinder
without Aspiration Metric / Us
Thread. chromium plated
Brass. - 1.8 CC & 2.2 CC



MSI-1743-14
3 Ring Cartridge Syringe for
Anesthesia - Non-Folding type
Available in Us / Metric Threads



MSI-1743-15
Articulated Syringe S-Type
Arpon Bayer for Anesthesia
Non folding. chromium
Plated Brass. Available in
Us / Metric Threads



MSI-1743-16
Aspirating Folding Cartridge
Syringe. chromium Plated
Brass. 1.8 CC & 2.2 CC
Available in Us/ Metric Thread



MSI-1743-17
Standard Syringe
1.8 ml
2.2 ml



MSI-1743-18

Aspirating Folding Syringe - strip Handle, single side folding type chromium Plated Brass. 1.8 CC & 2.2 CC Available in Us/ Metric Thread



MSI-1743-19

Pro Millex "A" (ASTRA) Type Aspirating Syringe with Metric /Us Threads. Non folding, retract by pulling thumb ring back. Insert cartridge, plunger end first into barrel. Backward motion will produce aspiration. 1.8 CC & 2.2 CC.



MSI-1743-20

Cork Screw Tip Aspirating syringe for anesthesia, single ring, Cylinder with Aspiration Metric / us Thread-chromium plated Brass. 1.8 CC & 2.2 CC.



MSI-1743-21

Cork Screw Tip Aspirating syringe for anesthesia; single ring, Cylinder with Aspiration Metric / us Thread-chromium plated Brass. 1.8 CC & 2.2 CC.



MSI-1743-22

Cartridge Syringe Cylinder with Aspiration, Metric /Us threads. Non folding- chromium Plated brass. 1.8 CC & 2.2 CC



Single ring aspirating syringe

MSI-1743-23

1.8 ml
2.2ml



MSI-1743-24


Cook-waite Type Self Aspirating Syringe, 1.8 CC & 2.2 CC for Anesthesia, Light weight Construction of Aluminum & stainless steel, available in Us / Metric threads & Sterilizable



MSI-1743-25

Self Aspirating Aspiject Syringe with saddle grip, for anesthesia stainless steel & High Tech polymer Construction which allows it to be completely Sterilizable

Ampullenspritze Cartridge Syringes



MSI-1743-26

Veterinary Metal Syringe made of copper & finished by chromium plating. High Quality that it will stand the test of time. Designed for one hand operation & Are Easy to Disassemble & clean. Disinfect it in boiled water. 10ml, 20ml, 30ml, 50 ml is Available.



MSI-1743-27

Detail: Component from (I to V) are made of copper & finished by chromium plating functions: Mainly Used for animal Disease Prevention & Treatment by Vet Not. The Disinfection in boiled water is applicable for syringes. High-pressure steam disinfection is prohibited.



MSI-1743-28

Veterinary Pistol Grip Pro-shot syringes Europlex Metal gun syringe / Pro-shot 50 cc is Ruggedly built from Highimpact plastic to be both Durable, Transparent, Graduated Barrel, which will not discolor. Simple dosage sleeve can be Easily Adjusted from 1 to 5 Ccs. Non-adjustable nylon plunger assembly & A Luer Lock Hub, extremely Accurate Pro-shot 5 cc Pistol-grip syringe combines outstanding features.



MSI-1743-29

Veterinary Syringe C-Type A Durable Plastic Syringe that will stand the test of time. Clear plastic Barrel, Black inscribed Graduations & orange Handles. Syringes have luer-lock luer-slip Needle Holder & are Autoclavable up to 150 Available in 10 ml, 20ml, 30ml, 50ml, & 100ml.



MSI-1743-29

Reiner Ear Syringe with 3 Cannulas & 1 Protecting Disk, Available in 50ml, 75ml, 100ml & 150ml.



MSI-1743-30

Aqua-fix Irrigating Self Filling Water Syringe (all Metal) for Ear.



MSI-1743-31

Ear Syringe with Spring Chrome plated

Chip Blowers Luftblasser



MSI-1744-00



MSI-1744-01



MSI-1745-00



MSI-1745-01



MSI-1746-01
without Valve



MSI-1746-02
with Valve

