

Gusslegierung Flaschen & Presse Denture Alloy Flasks & Press



MSI-2226-00

Upper / Lower Denture Flask & Bench Press
Denture Flask Set featur forged Brass alloy
Construction. Fitted with Hardened stainless steel
alloy Rear Tapered Guide pins. Flask Dimensions:
Long x 3^{3/4}Wide x 2^{1/4}High
(10cm x 10cm x 6cm) Unit wweight: 3.850kg including
(upper / Lower Flask & Bench press)



MSI-2227-00

Small Denture Flask & Bench Press
Denture Flask Set featur forged Brass alloy
Construction. Fitted with Hardened stainless steel
alloy Rear Tapered Guide pins. Flask Dimensions:
4"Long x 3^{3/4}Wide x 2^{1/4}High
(10cm x 10cm x 6cm) Unit wweight: 2.425kg including
(Flask & Bench press)



MSI-2228-00

Flat Bench Press for 2 Denture
Flasks (upper + lower).
Dimension: 20 x 12.7 x 8cm &
Unit weiht is 2.17kg



MSI-2229-00

Alloy Brass Round Bench Press
for 2 Flask. Dimension: 7"x5"x 3^{1/2}"
& Unit Weight 1.900kg



MSI-2230-00

Small Denture Bench Press
for single Flask. Brass
Alloy Construction.
Dimension: 5"x 3"x 4^{3/4}"
& Unit Weight 1.500kg



MSI-2231-00

Large Denture Flask &
Bench Press. Alloy Construction
Dimension: Flask Dimension:
5"x 3^{1/2}"x 2^{1/2}"
& Weight 1.100kg Each
Unit weight 5.750kg
Icluding (3 Flasks & Bench Press)



MSI-2232-00

Large Bench Press for
3 Flasks. Featured forged
Brass Alloy construction.
Dimension: 29 x 14 x 9cm
Unit weight: weight 2.800kg

Gusslegierung Flaschen & Presse Denture Alloy Flasks & Press



MSI-2233-00
Round Denture Flask
are made in superior
Gunmetal & Brass
The round Flasks are Appropriate
Dimension: 4" x 4 1/2" x 2 1/2" and
Weight 1.100kg



MSI-2234-00
Crown Flask for Dental Tooth,
Forged Brass Alloy construction
Unit weight: 365 gram



MSI-2235-00
One Side Flat Brass Flask,
New Shape forged brass
Alloy construction,
Dimension: 2 1/4" x 3 1/2" x 4 3/4"
& Unit weight 955 gram



MSI-2236-00
Mortimer Silverline Alloy Flask,
Forged Silver Alloy construction,
Dimension: 4 1/4" Long x 4"
Wide x 2 3/4" High (11 x 10 x 7cm)
Unit weight = 750 gram



MSI-2237-00
Forged Brass Toothed
Dental Alloy Flask for
Single tooth,
Dimension: 5.0 x 5.0 x
4.0cm &
Unit weight 270gram



MSI-2238-00
Square Dental Alloy Flask,
Pressure Die Casted, made of brass
Dimension: 7.2 x 8.4 x 5.0 cm
Unit weight 1.300kg



MSI-2239-00
Upper Denture Flask forged
Brass Alloy construction Fitted
with Hardened stainless steel
Alloy rear Tapered guide pin:
Dimension: 4" Long x 3 3/4"
wide x 2 1/4" High (10cm x 10cm x 6cm)
& weight: 0.920kg



MSI-2240-00
Lower Denture Flask forged
Brass Alloy construction. Fitted
with Hardened stainless steel
Alloy rear Tapered guide pins.
Dimension: 4" Long x 3 3/4"
wide x 2 1/4" High (10cm x 10cm x 6cm)
& weight: 0.920kg

Gusslegierung Flaschen & Presse
Denture Alloy Flasks & Press



MSI-2241-00

Upper Denture Flask forged Silver Alloy construction. Fitted with Hardened stainless steel Alloy rear Tapered guide pins. Dimension: 4" Long x 3 3/4" wide x 2 1/4" High (10cm x 10cm x 6cm)



MSI-2242-00

Lower Denture Flask forged Silver Alloy construction. Fitted with Hardened stainless steel Alloy rear Tapered guide pins. Dimension: 4" Long x 3 3/4" wide x 2 1/4" High (10cm x 10cm x 6cm)




MSI-2243-00

Upper / Lower Hanau Denture Flask, forged Brass Alloy Construction. Fitted with Hardened stainless steel Alloy rear Tapered guide pins. Dimension: 10 x 8 x 16 cm wide x 2 1/4" Long x 3 3/4"



MSI-2244-00

Silverline Alloy Flask, forged Silver Alloy construction, Dimension: 4 1/4" Long x 4" wide x 2 3/4" High (11 x 10 x 7cm) Unit weight = 760 gram



MSI-2245-00

Butter Dish Silver Alloy Denture Flask. Dimension: 10 x 8 x 16 cm Unit weight = 760 gram



MSI-2246-00

Flask Knock-out Disc Forged Brass Alloy construction



MSI-2247-00
Gic Applicator



MSI-2248-00
Gic Activator

**Maisschneider
Corn Cutter**



MSI-2249-01
Corn Cutter



MSI-2249-02
Corn Cutter - Plastic Handle



MSI-2249-03
Corn Cutter - Plastic Handle



MSI-2249-04
Corn Cutter - Plastic Handle



MSI-2250-00
Universal Bur
Changer Wrench



MSI-2251-00
Dentaire Bur
Changer wrench



MSI-2252-00
Pocket Instruments
Holder



MSI-2253-00
Bow Caliper 8
(compasse)



MSI-2254-00
Diamond Disc Guard
with adapter



MSI-2255-00
Caliper Eight Gauge



MSI-2256-00
Optical Inspection stand



MSI-2257-00
Matrix Retainer Elliott
Separator



MSI-2258-00
Gem Holder
Bead Holder



MSI-2259-00
Tongue Depressor
Small, Large



MSI-2260-00
Ring Bending Pliers
Swiss Model 190 mm



MSI-2261-00
Mesh Bur Holder



MSI-2262-00
Stripholder Clips



MSI-2263-00
X-Ray Film
Holder



MSI-2264-00
Radiography Clip



MSI-2265-00
X-Ray Film hangers Multiple Clips
Available in 6, 8, 10, 12, & 14 Clips



MSI-2266-00
Feline Mouth Gag for
cats with Nylon Inserts
Available in 2 sizes
Small 2.5" Large 3"



MSI-2267-00
Canine Spring Action
Mouth Gag with Nylon Inserts
Available in 3 Sizes Small 4.5"
Medium 5.28" - Large 7"



MSI-2268-01
Mouth Mirror Keychains



MSI-2268-02
Mouth Mirror Keychains



MSI-2268-03
Mouth Mirror Keychains



MSI-2268-04
Mouth Mirror Keychains



MSI-2268-05
Tweezer Keychains



MSI-2268-06
Dental Tooth Keychains



MSI-2268-07
Extracting Forceps Keychains



MSI-2268-08
Dental Caliper Keychains



MSI-2268-09



MSI-2268-10
Dental Syringe
Keychains



MSI-2268-11
Dental Syringe
Keychains



MSI-2268-12
Impression Tray Keychains



MSI-2269-01
Dental Forceps Set
Childs



MSI-2269-04
Dental Forceps Set



MSI-2269-02
Rubber Dam Set



MSI-2269-05



MSI-2269-06



MSI-2269-03
Rubber Dam Set



MSI-2269-07
Root Elevator Set

Luxator Satz Von 8 stück
Luxator Set of 8 Pcs



MSI-2270-01



MSI-2270-02



MSI-2270-03



MSI-2270-04



MSI-2270-05



MSI-2270-06



MSI-2270-07



MSI-2270-08