


MACINTOSH FIBER OPTIC



MSI-8729-29



ADVANTAGES

Fiber optic blades with German Bulb
Fiber optic illumination
assembly is detachable
and is easy to clean.
Maintain and re-polished.
Thousands of fiber strings
are firmly secured in stainless steel
pipe cover which saves
from damage or breakage.

MILLER FIBER OPTIC



MSI-8730-30



MILLER CONVENTIONAL



MSI-8731-01



Fig.4



Fig.3



Fig.2



Fig.1

ADVANTAGES

Handles

Brass
Universal Fitting with
Non Magnetic pin

Blades

Stainless Steel
Non Magnetic
with german bulb

MACINTOSH CONVENTIONAL



MSI-8732-02



Fig.4



Fig.3



Fig.2



Fig.1



DIAGNOSTIC ENT SET



MSI-8733-33



OTOSCOPE SET WITH CASE



MSI-8734-34

Laryngoscope Set, Handle & Blades



MSI-8735-00



Miller Baby
MSI-8736-00

FIBEROPTIC HANDLES



Fig.1



Fig. 2



Fig. 3



Fig. 4

MSI-8738-00

ADVANTAGES

Choice of Xenon & Halogen lamps available.
All types of standard lamp threads available.
Easy to replace lamp
Handles fitting is compatible

ADVANTAGES

High quality long life handles
Handles fitting is compatible



MSI-8737-00
Magil Set of 4 Blade
Size : 1,2,3,4 1Medium Handle
1 Small & 1Large spare Lamp packed in
Strong Plastic Casing with Foam Inserting

CONVENTIONAL HANDLES



Fig.1



Fig. 2



Fig. 3



Fig. 4

MSI-8738-39

TONGUE DEPRESSOR



MSI-8738-A

MINI OTOSCOPE FIBER OPTICS

PEN TORCH



MSI-8738-B



MSI-8738-C



MSI-8738-D



MSI-8738-E

Available in colour



MSI-8738-F

MSI-8738-G

MSI-8738-H

METAL BODY OTOSCOPIES



MSI-8738-J

MSI-8738-K




MSI-8738-L

OPHTHALMOSCOPIES




MSI-8738-M

MSI-8738-N




| Ømm | Fergusson |
|-------|----------------|
| 2 | 6 MSI-8739-06 |
| 2 1/2 | 7 MSI-8739-07 |
| 2 2/3 | 8 MSI-8739-08 |
| 3 | 9 MSI-8739-09 |
| 3 1/2 | 10 MSI-8739-10 |
| 4 | 11 MSI-8739-12 |
| 4 3/4 | 12 MSI-8739-14 |

17cm/6 3/4"



| Ømm | Frazier |
|-------|----------------|
| 2 | 6 MSI-8740-06 |
| 2 1/2 | 7 MSI-8740-07 |
| 2 2/3 | 8 MSI-8740-08 |
| 3 | 9 MSI-8740-09 |
| 3 1/2 | 10 MSI-8740-10 |
| 4 | 12 MSI-8740-12 |
| 4 3/4 | 14 MSI-8740-14 |
| 5 1/4 | 16 MSI-8740-16 |

17cm/6 3/4"



Yasargil
Luer boring
Luer langweilig
ennuyeux Luer
aburrido Luer

| usable length cm | Ømm diam |
|------------------|-----------------|
| 13 | 1.5 MSI-8741-13 |
| 15 | 1.5 MSI-8741-15 |
| 18 | 1.5 MSI-8741-18 |
| 22 | 1.5 MSI-8741-22 |
| 13 | 2.0 MSI-8742-13 |
| 15 | 2.0 MSI-8742-15 |
| 18 | 2.0 MSI-8742-18 |
| 22 | 2.0 MSI-8742-22 |
| 13 | 2.5 MSI-8743-13 |
| 15 | 2.5 MSI-8743-15 |
| 18 | 2.5 MSI-8743-18 |
| 22 | 2.5 MSI-8743-22 |
| 13 | 3.5 MSI-8744-13 |
| 15 | 3.5 MSI-8744-15 |
| 18 | 3.5 MSI-8744-18 |
| 22 | 3.5 MSI-8744-22 |
| 18 | 4.5 MSI-8744-18 |

| Ømm | |
|-----|-------------|
| 1.7 | MSI-8745-17 |
| 2.5 | MSI-8746-25 |
| 3.5 | MSI-8747-35 |

malleable
formbar
malleable
maleable

Yasargil

Anaesthesia Instruments

Suction Tubes
Saugrohre

Tubes d'aspiration
tubos de succión



malleable, black
Temperguss, schwarz
malleable, noire
maleable, negro

∅

| | |
|------|-------------|
| 2 mm | MSI-8748-01 |
| 3 mm | MSI-8748-02 |
| 4 mm | MSI-8748-03 |
| 5 mm | MSI-8748-04 |

19.5 cm/7"

Brain Suction Tubes
Gehirn Saugschläuche
Tubes d'aspiration du cerveau
Los tubos de succión del cerebro



∅

Plester

| | |
|--------|-------------|
| 1.5 mm | MSI-8750-15 |
| 2.0 mm | MSI-8750-20 |
| 2.5 mm | MSI-8750-25 |
| 3.0 mm | MSI-8750-30 |

19.5cm/7 3/4"



Adson
∅4 mm

Adson
MSI-8749-50
21 cm/8 1/4"



∅

Frazier

| | |
|--------|-------------|
| 2.0 mm | MSI-8751-20 |
| 2.5 mm | MSI-8751-25 |
| 3.0 mm | MSI-8751-30 |

13cm/5"

Ø 8 mm



Poole
MSI-8752-59
22cm/8¾"

Ø 10 mm



Poole
MSI-8753-60
22cm/8¾"


Ø 6 mm



Baby-Poole
MSI-8754-61
20cm/8"



Finsterer
MSI-8755-62
30cm/12"



17.5 cm/7"
MSI-8756-63



Fig. 3
MSI-8757-03



Fig. 2
MSI-8757-02



Fig. 1
MSI-8757-01



Fig. 0
MSI-8757-00



Yankauer
MSI-8758-53



Yankauer
MSI-8759-54
27cm/10½"

with stopcock
mit Hahn
avec robinet d'arrêt
con llave de paso



Yankauer
MSI-8760-55
27cm/10½"




Andrew-Pynchon
MSI-8761-56
23 cm/9"



Nuboer
MSI-8762-57
21cm/8¼"




Nuboer
MSI-8763-58
21cm/8¼"



Magill

MSI-8764-01=20cm/8"
MSI-8764-02=20cm/8"
MSI-8764-03=20cm/8"



**Yankauer universal
Complete Set
MSI-8765-65**



**Yankauer
MSI-8766-34
34.5cm/13 1/4"**



**Simpson-Smith
MSI-8767-67**





**Millin
MSI-8768-30
30cm/12"**

Anaesthesia Instruments

Nasal and Ear Cannulas
Nasen-und Ohr-Kanülen

Nasen-und Ohr-Kanülen...
Cánulas nasales y del oído






Ayre
 MSI-8775-01 Small Small Medium Large
 MSI-8775-02 Medium klein Medium groß
 MSI-8775-03 Large petit moyen grand
 pequeño medio grande




Cobb
 MSI-8776-01 Small Small Medium Large
 MSI-8776-02 Medium klein Medium groß
 MSI-8776-03 Large petit moyen grand
 pequeño medio grande



Mitchell
 MSI-8777-18
 MSI-8777-20
 38mm



Hamilton Bailey
 MSI-8778-01 Infant Infant Child Adult
 MSI-8778-02 Child Säugling Kind Erwachsene
 MSI-8778-03 Adult nourrisson enfant adulte
 infantil niño adulto



Magill
 Connection Set
 MSI-8779-00



| Magill (Nasal) | | Magill (Nasal) | |
|----------------|--------|----------------|--------|
| | mm dia | | mm dia |
| MSI-8780-01 | 5.5 | MSI-8781-01 | 5.5 |
| MSI-8780-02 | 6.5 | MSI-8781-02 | 6.5 |
| MSI-8780-03 | 7.5 | MSI-8781-03 | 7.5 |
| MSI-8780-04 | 8.0 | MSI-8781-04 | 8.0 |
| MSI-8780-05 | 8.5 | MSI-8781-05 | 8.5 |
| MSI-8780-06 | 9.5 | MSI-8781-06 | 9.5 |
| MSI-8780-07 | 10.0 | MSI-8781-07 | 10.0 |
| MSI-8780-08 | 10.5 | MSI-8781-08 | 10.5 |