


Scheme of Dimensions for Specula



	a mm	b mm
MSI-9800-01	75	32
MSI-9800-02	85	35
MSI-9800-03	100	37

Cusco



	a mm	b mm
MSI-9801-01	75	32
MSI-9801-02	85	35
MSI-9801-03	100	37

Cusco



For Virgins	a mm	b mm
MSI-9802-00	75	17

Cusco



For Virgins	a mm	b mm
MSI-9803-00	75	17

Cusco



	a mm	b mm
MSI-9804-01	90	23-25
MSI-9804-02	100	25-27
MSI-9804-03	110	27-30

Cusco
Swiss Pattern

Gynecology Instruments

Vaginal Speculas
Vaginal Speculas

Speculas vaginales
especulos vaginales




Collin-Landau

	a mm	b mm
MSI-9805-01	85	30
MSI-9805-02	100	35
MSI-9805-03	110	40



Collin

	a mm	b mm
MSI-9806-01	85	30
MSI-9806-02	100	35
MSI-9806-03	110	40



Trelat

	a mm	b mm
MSI-9807-01	85	33
MSI-9807-02	95	35
MSI-9807-03	115	43

Gynecology Instruments

Vaginal Speculas For Laser Surgery
Vaginal Speculas Für Laser-Chirurgie

Speculas vaginaux pour la chirurgie laser
Especulos vaginales para la cirugía láser

MEDSIGN INSTRUMENTS



Grave (Standard)

	a mm	b mm
MSI-9808-01	75	20
MSI-9808-02	95	30
MSI-9808-03	115	35
MSI-9808-04	125	35

Black




Grave	a mm	b mm
MSI-9809-01	75	20
MSI-9809-02	95	30
MSI-9809-03	115	35



Insulated

Grave	a mm	b mm
Angled 135°	101	35
MSI-9810-00		



Grave

Right Side open
Right Side offen
Côté droit ouvert
Lado derecho abierto

Grave

Right Side open Angled
Right Side öffnen Angled
Côté droit à angle ouvert
Lado derecho abrir en ángulo

	a mm	b mm
MSI-9811-01	101	35
MSI-9811-02		

Gynecology Instruments

Vaginal Speculas
Vaginal Speculas

Speculas vaginales
especulos vaginales

Lateral Specula



MSI-9812-00



Pederson	a mm	b mm
MSI-9813-01	75	13
MSI-9813-02	100	22
MSI-9813-03	120	25



Sayfert	a mm	b mm
MSI-9814-01	100	17
MSI-9814-02	1010	33



MSI-9815-01	small, klein, petit, pequeño
MSI-9815-02	medium, medeum, medeum
MSI-9815-03	large, groß, grand, grande

Scheme of Dimensions for Specula
Schema Abmessungen für Spekula
Schéma des dimensions de spéculums
Esquema de Medidas para espejuelos



a
Length of Blades
ein Länge der Klingen
un Longueur des lames
un Duración de Blades



b
Width of Blades
b Breite der Klingen
b Largeur des lames
b Ancho de Espadas



Sims

MSI-9816-01 63 x 32 mm
MSI-9816-02 76 x 32 mm
MSI-9816-03 89 x 38 mm



Sims

	Fig	a mm	b mm
MSI-9817-00			
MSI-9817-01	1	60	25+30
MSI-9817-02	2	70	30+35
MSI-9817-03	3	80	35+40



Jackson

Fig.1 75 x 35 MSI-9818-01
Fig. 2 90 x 35 MSI-9818-02
Fig. 3 100 x 35 MSI-9818-03

	Fig.	a mm	b mm
MSI-9818-04	4	75 X 40	
MSI-9818-05	5	90 X 40	
MSI-9818-06	6	100 X 40	



Eastman

	Fig.	a mm	b mm
MSI-9819-01	1	76	38
MSI-9819-02	2	89	38
MSI-9819-03	3	101	38



For children
für Kinder
pour les enfants
Para los niños

Set
setzen
fixer
establecer

Kristeller
MSI-9820-00
70x15 mm



Specula

Retractor

Kristeller
(Standard)

Set of 1 specula and 1 retractor
Set 1 Spekula und 1 Retraktor
Set de 1 spéculum et 1 enrouleur
1 Juego de espéculos y retractor 1

Fig. 1 MSI-9821-00

Fig. 2 MSI-9822-00

Fig. 3 MSI-9823-00

Specula

	Fig.	a mm	b mm
MSI-9821-01	1	70	27
MSI-9822-01	2	80	30
MSI-9823-01	3	90	36

Retractors

	Fig.	a mm	b mm
MSI-9821-02	1	75	23
MSI-9822-02	2	85	26
MSI-9823-02	3	100	32



With conically shaped blades for virgins
Mit konisch geformten Klingen für Jungfrauen
Avec de forme conique lames pour vierges
Con forma cónica cuchillas para vírgenes

Kristeller



With extra long blades
Mit extra langen Klingen
Avec de longues lames supplémentaires
Con extra largas hojas

**Kristeller
(Standard)**

Set of 1 specula and 1 retractor
Set 1 Spekula und 1 Retraktor
Set de 1 spéculums et 1 enrouleur
1 Juego de especúlos y retractor 1

Fig. 00 MSI-9824-00

Fig. 0 MSI-9825-00

Specula

	Fig.	a	b
		mm	mm
MSI-9824-01	00	110	18
MSI-9825-01	0	110	21

Retractors

	Fig.	a	b
		mm	mm
MSI-9824-02	00	115	14
MSI-9825-02	0	115	17

Set of 1 specula and 1 retractor
Set 1 Spekula und 1 Retraktor
Set de 1 spéculums et 1 enrouleur
1 Juego de especúlos y retractor 1

Fig. 1 MSI-9826-00

Fig. 2 MSI-9827-00

Fig. 3 MSI-9828-00

Specula

	Fig.	a	b
		mm	mm
MSI-9826-01	1	110	27
MSI-9827-01	2	110	30
MSI-9828-01	3	110	36

Retractors

	Fig.	a	b
		mm	mm
MSI-9826-02	1	115	23
MSI-9827-02	2	115	26
MSI-9828-02	3	115	32

Gynecology Instruments

Vaginal Speculas
Vaginal Speculas

Speculas vaginales
espéculos vaginales



With flat hollow specula
Mit flachen Talkessel-Spekula
Avec spéculums à plat creux
Con espéculos cuenca llana

Set of 1 specula and 1 retractor
Set 1 Spekula und 1 Retraktor
Set de 1 spéculums et 1 enrouleur
1.Juego de espéculos y retractor 1

Fig. 1 MSI-9829-00

Fig. 2 MSI-9830-00

Fig. 3 MSI-9831-00

Specula

	Fig.	a	b
		mm	mm
MSI-9829-01	1	105	26
MSI-9830-01	2	105	28
MSI-9831-01	3	105	33

Retractors

	Fig.	a	b
		mm	mm
MSI-9829-02	1	115	22
MSI-9830-02	2	115	24
MSI-9831-02	3	115	29



	a	b
	mm	mm
MSI-9832-20	100	20
MSI-9832-25	100	25
MSI-9832-30	100	30
MSI-9832-35	100	35
MSI-9832-40	100	40
MSI-9832-20	130	20
MSI-9832-25	130	25
MSI-9832-30	130	30
MSI-9832-35	130	35
MSI-9832-40	130	40

Matzenauer



Braun
MSI-9833-18
18 cm/7"
56x13 mm




Braun
MSI-9834-18
18 cm/7"
60x10 mm



Ulrich
MSI-9835-19
19.5 cm/7 1/2"
80x20 mm



MSI-9836-19
19 cm/7 1/2"
85x17 mm



Kallmorgen

Set

MSI-9837-00	Fig. 1
MSI-9838-00	Fig. 2

Specula

	Fig.	a	b
		mm	mm
MSI-9837-01	1	70	40
MSI-9838-01	2	90	40


Retractors

	Fig.	a	b
		mm	mm
MSI-9837-02	1	70	33
MSI-9838-02	2	90	33

Gynecology Instruments


Vaginal Speculas
Vaginal Speculas

Speculas vaginales
especulos vaginales



Mathieu

Specula	Fig.	a	b
		mm	mm
MSI-9839-01	1	100	13/23
MSI-9840-01	2	100	16/30



Mathieu

Retractors	Fig.	a	b
		mm	mm
MSI-9839-02	1	100	11/21
MSI-9840-02	2	100	14/27

Set	Fig.	MSI-9839-00
		MSI-9840-00
	1	MSI-9839-00
	2	MSI-9840-00




Doyen

Fig.	a	b
	mm	mm
MSI-9841-01	75	45
MSI-9841-02	90	45
MSI-9841-03	110	45



Doyen
24 cm/9 1/2"
slightly-concave
leicht konkav-
légèrement concave
ligeramente cóncava

	Fig.	a mm	b mm
MSI-9843-01	1	55	35
MSI-9843-02	2	85	35
MSI-9843-03	3	115	35
MSI-9844-01	1	55	45
MSI-9844-02	2	85	45
MSI-9844-03	3	115	45
MSI-9844-04	4	160	45
MSI-9844-05	5	180	45
MSI-9845-01	1	55	60
MSI-9845-02	2	85	60
MSI-9845-03	3	115	60
MSI-9845-04	4	160	60
MSI-9845-05	5	180	60



Doyen
24 cm/9 1/2"
Medium-concave
Medium-konkaven
Moyennes concave
Medio cóncavo

MSI-9846-01
MSI-9846-02
MSI-9846-03
MSI-9847-01
MSI-9847-02
MSI-9847-03
MSI-9848-01
MSI-9848-02
MSI-9848-03



Doyen
25 cm/10"

	a mm	b mm
MSI-9849-01	60	30
MSI-9849-02	70	30
MSI-9849-03	80	30



Mangiagalli
MSI-9850-00
40 x 80 mm

Gynecology Instruments

Vaginal Speculas
Vaginal Speculas

Speculas vaginales
especulos vaginales



Auvard

with removable weight
mit abnehmbarem Gewicht
avec un poids amovible
con peso extraible

MSI-9851-01
80x38mm
MSI-9851-02
105x42mm

Auvard

with fixed weight
mit festem Gewicht
avec un poids fixe
con peso fijo

MSI-9852-01
80x38 mm
MSI-9852-02
105x42 mm



Landau

29 cm/11½"

	a	b
	mm	mm
MSI-9853-01	100	28
MSI-9853-02	100	38



Breisky

	a	b
	mm	mm
MSI-9854-20	100	20
MSI-9854-25	100	25
MSI-9854-30	100	30
MSI-9854-35	100	35
MSI-9854-40	100	40
MSI-9855-20	130	20
MSI-9855-25	130	25
MSI-9855-30	130	30
MSI-9855-35	130	35
MSI-9855-40	130	40
MSI-9856-20	160	20
MSI-9856-25	160	25
MSI-9856-30	160	30
MSI-9856-35	160	35
MSI-9856-40	160	40
MSI-9857-20	180	20
MSI-9857-25	180	25
MSI-9857-30	180	30
MSI-9857-35	180	35
MSI-9857-40	180	40

Gynecology Instruments

Vaginal Speculas and Retractors Speculas vaginaux et rétracteurs
 Vaginal Speculas und Retraktoren Espéculos vaginales y retractoros

MEDSIGN INSTRUMENTS



Bozemann
MSI-9858-00

Complete



MSI-9858-05
80 x 30 mm



MSI-9858-01
65 x 30 mm



MSI-9858-06
85 x 33 mm



MSI-9858-02
75 x 35 mm



MSI-9858-07
100 x 38 mm



MSI-9858-03
85 x 35 mm



MSI-9858-08
110 x 42 mm



MSI-9858-04
95 x 35 mm



O'Sullivan-O'Connor
MSI-9859-00

Complete
vervolständigen
completter
completar



	Fig.	a mm	b mm
MSI-9860-01	1	45	17
MSI-9860-02	2	90	37
MSI-9860-03	3	90	45

Gynecology Instruments

Uterine Dilators
Uterus-Dilatatoren

Dilatateurs utérins
dilatadores uterinos



Hegar

mm Ø	mm Ø
1.0 MSI-9861-01	1.5 MSI-9862-01
2.0 MSI-9861-02	2.5 MSI-9862-02
3.0 MSI-9861-03	3.5 MSI-9862-03
4.0 MSI-9861-04	4.5 MSI-9862-04
5.0 MSI-9861-05	5.5 MSI-9862-05
6.0 MSI-9861-06	6.5 MSI-9862-06
7.0 MSI-9861-07	7.5 MSI-9862-07
8.0 MSI-9861-08	8.5 MSI-9862-08
9.0 MSI-9861-09	9.5 MSI-9862-09
10.0 MSI-9861-10	10.5 MSI-9862-10
11.0 MSI-9861-11	11.5 MSI-9862-11
12.0 MSI-9861-12	12.5 MSI-9862-12
13.0 MSI-9861-13	13.5 MSI-9862-13
14.0 MSI-9861-14	14.5 MSI-9862-14
15.0 MSI-9861-15	15.5 MSI-9862-15
16.0 MSI-9861-16	16.5 MSI-9862-16
17.0 MSI-9861-17	17.5 MSI-9862-17
18.0 MSI-9861-18	18.5 MSI-9862-18
19.0 MSI-9861-19	19.5 MSI-9862-19
20.0 MSI-9861-20	
21.0 MSI-9861-21	
22.0 MSI-9861-22	
23.0 MSI-9861-23	
24.0 MSI-9861-24	
25.0 MSI-9861-25	
26.0 MSI-9861-26	
27.0 MSI-9861-27	
28.0 MSI-9861-28	
29.0 MSI-9861-29	
30.0 MSI-9861-30	

Pieces per set	set
MSI-9861-40 14	4.0-17.0 mm In canvas
MSI-9861-50 14	4.0-17.0 mm In metal case
MSI-9861-60 26	1.0-26.0 mm in canvas
MSI-9861-70 26	1.0-26.0 mm in metal case



Hegar


Ø mm	
1.0/2.0 MSI-9863-02	
3.0/4.0 MSI-9863-04	
5.0/6.0 MSI-9863-06	
7.0/8.0 MSI-9863-08	
9.0/10.0 MSI-9863-10	
11.0/12.0 MSI-9863-12	
13.0/14.0 MSI-9863-14	
15.0/16.0 MSI-9863-16	
17.0/18.0 MSI-9863-18	
19.0/20.0 MSI-9863-20	
21.0/22.0 MSI-9863-22	
23.0/24.0 MSI-9863-24	
25.0/26.0 MSI-9863-26	

Pieces per set	set
MSI-9863-40 08	3.0/4.0-17.0/18.0 mm
MSI-9864-60 13	1.0/2.0-25.0/26.0 mm
MSI-9863-50 08	3.0/4.0-17.0/18.0 mm
MSI-9864-70 13	1.0/2.0-25.0/26.0 mm



Pratt

Fig.	mm	
13/15	4.5/5.0	MSI-9865-15
17/19	5.6/6.3	MSI-9865-19
21/23	7.0/7.6	MSI-9865-23
25/27	8.3/9.0	MSI-9865-27
29/31	9.6/10.3	MSI-9865-31
33/35	11.0/11.6	MSI-9865-35
37/39	12.3/13.0	MSI-9865-39
41/43	13.6/14.3	MSI-9865-43



Set
MSI-9865-00



Hank
U20-9866-00

Fig.	mm	
9/10	4.5/5	MSI-9866-10
11/12	5.5/6	MSI-9866-12
13/14	6.5/7	MSI-9866-14
15/16	7.5/8	MSI-9866-16
17/18	8.5/9	MSI-9866-18
19/20	9.5/10	MSI-9866-20



Set

Gynecology Instruments

Uterine Tenaculum Forceps, Amniotome Tenaculum utérine Forceps, amniotomie
 Uterine Hákchen Pinzette, Amniotome Uterina tenáculo pinzas, Amniotome



Sims
 MSI-9867-33
 33 cm/13"



Sims
 MSI-9868-33
 33 cm/13"



Walleix
 MSI-9869-29
 29.5 cm/11½"



Braun
 MSI-9870-27
 27.5 cm/11¾"

Gynecology Instruments

Uterine Sounds, Uterine Depressor
Uterus-Sounds, Uterus-Depressor

Sons utérins, de l'utérus dépresseurs
Los sonidos uterinos, útero depresores



Schröder
MSI-9871-25
25.5 cm/10 "



Duply
MSI-9872-25
25.5 cm/10 "



Braun
MSI-9873-25
25 cm/10 1/4"



Beacham
MSI-9874-26
26 cm/10 1/4"

Gynecology Instruments

Uterine Tenaculum Forceps
Uterine Tenakulum Zange

Utérine Forceps tenaculum
Pinza uterina tenáculo



Schroder



MSI-9875-25
25 cm/9 3/4"



MSI-9876-23
23.5 cm/9 1/4"



Skene
MSI-9877-24
24 cm/9 1/2"



Teale
MSI-9878-23
23.5 cm/9 1/4"

Gynecology Instruments

Uterine Tenaculum Forceps
Uterine Tenakulum Zange

Utérine Forceps tenaculum
Pinza uterina tenáculo

MEDSIGN INSTRUMENTS



Iowa



MSI-9879-26
26 cm/10 1/4"



Jacobs



MSI-9883-21
21 cm/8 1/4"



Museux



6 mm
MSI-9880-24
24 cm/9 1/2"



6 mm
MSI-9881-24
24 cm/9 1/2"



6 mm
MSI-9882-24
24 cm/9 1/2"

Gynecology Instruments

Uterine Tenaculum Forceps
Uterine Tenakulum Zange

Utérine Forceps tenaculum
Pinza uterina tenáculo



MSI-9884-26
26.5 cm/10 1/2"



MSI-9885-16
16 cm/10 1/4"

Braun



MSI-9886-21
21 cm/8 1/4"



Czerny

Gynecology Instruments

Uterine Polypus Forceps
Uterine Polypus Zange

Utérine Polypus Forceps
Pinza uterina Polypus



MSI-9887-32
32 cm/12 1/4"



MSI-9888-27
27 cm/10 1/4"



MSI-9889-25
25 cm/9 3/4"




Gynecology Instruments

Uterine Seizing Forceps
Uterine ergriff Zange

Utérins Forceps Saisir
Uterinos Pinzas Aprovechar



Dartigues
MSI-9890-25
28 cm/11"



Somer
MSI-9891-25
22.5 cm/8 3/4"



Collin
MSI-9892-26
26 cm/10 1/4"

Gynecology Instruments

Myoma Screws, Myoma Knife
Schrauben Myom, Myome Knife

Vis myome, couteaux myomes
Tomillos, mioma Mioma Knife




Doyen
MSI-9893-17
17 cm/6 3/4"



Doyen
MSI-9894-14
14.5 cm/5 3/4"



Segand



MSI-9895-28
28 cm/11"

Gynecology Instruments

Suction Curettes, Uterine Biopsy Curettes Curettes d'aspiration, de l'utérus biopsie Curettes
 Saug-Küretten, Uterus-Biopsie Küretten Curetas de aspiración, biopsia uterina Curetas



Novak



MSI-9896-24
 25.5cm/8 3/4"
 4.0 mm ø



Novorkian



MSI-9897-30
 30cm/11 3/4"
 6.0 mm ø



Munchen



MSI-9898-27
 27 cm/10 3/4"
 4.8 mm ø



MSI-9899-23
23.5 cm/9 1/4"

Ayre




MSI-9900-23
22 cm/9 3/4"

Emmett



MSI-9901-21
21.5cm/8 1/2"

Simon



MSI-9902-21
21.5cm/8 1/2"

Simon

Gynecology Instruments

Uterine Biopsy Forceps
Uterine Biopsiezangen

Pince à biopsie utérine
Pinzas de biopsia uterina



Schubert



MSI-9903-21
27.5 cm/8 1/2"



MSI-9903-28
28 cm/11"



MSI-9904-21
21 cm/8 1/4"



Faure



MSI-9905-24
24 cm/8 1/2"



Gellhorn

Gynecology Instruments

Uterine Biopsy Forceps
Uterine Biopsiezangen

Pince à biopsie utérine
Pinzas de biopsia uterina

MEDSIGN INSTRUMENTS



MSI-9906-21
21 cm/8 1/4"



MSI-9906-23
23.5 cm/9 1/4"



MSI-9908-27
27 cm/10 3/4"




Thoms-Gaylor



Schubert



Van Doren




MSI-9907-25
25.5 cm/10"

Gynecology Instruments

Uterine Biopsy Forceps
Uterine Biopsiezangen

Pince à biopsie utérine
Pinzas de biopsia uterina



Gynecology Instruments

Uterine Biopsy Forceps
Uterine Biopsiezangen

Pince à biopsie utérine
Pinzas de biopsia uterina



MSI-9911-25
25.5 cm/10"

Van Doren












MSI-9912-28
28 cm/11"

Alexander

Gynecology Instruments

Uterine Biopsy Forceps, IUP Grasping Forceps Utérine pinces à biopsie, pinces à préhension IUP
 Uterine Biopsiezangen, IUP Faßzange Uterina pinzas de biopsia, Pinzas IUP Agarre



		
	9,6 x 6,8 mm MSI-9915-01 24,5 cm/9 3/4"	9,6 x 6,8 mm MSI-9916-01 30 cm/11 3/4"
		
	11,6 x 8,3 mm MSI-9915-02 24,5 cm/9 3/4"	11,6 x 8,3 mm MSI-9916-02 30 cm/11 3/4"
		
	13,2 x 10,0 mm MSI-9915-03 24,5 cm/9 3/4"	13,2 x 10,0 mm MSI-9916-03 31 cm/12 1/4"
		
	16,0 x 11,6 mm MSI-9915-04 24,5 cm/9 3/4"	16,0 x 11,6 mm MSI-9916-04 31cm/12 1/4"
		
	18,0 x 12,6 mm MSI-9915-05 25 cm/9 3/4"	18,0 x 12,6 mm MSI-9916-05 31 cm/12 1/4"
		
	20,0 x 13,8 mm MSI-9915-06 25 cm/9 3/4"	20,0 x 13,8 mm MSI-9916-06 31 cm/12 1/4"
		

Gynecology Instruments

Uterine Curettes Curettes utérines
 Uterus-Küretten Curetas uterinas



Recamier
 MSI-9917-05
 MSI-9917-40



Recamier
 MSI-9918-05
 MSI-9918-40



Sims
 MSI-9919-05
 MSI-9919-40



Sims
 MSI-9920-05
 MSI-9920-40

Fig	mm	MSI-9917-05	MSI-9918-05	MSI-9919-05	MSI-9920-05
00	5	MSI-9917-05	MSI-9918-05	MSI-9919-05	MSI-9920-05
0	6	MSI-9917-06	MSI-9918-06	MSI-9919-06	MSI-9920-06
1	7	MSI-9917-07	MSI-9918-07	MSI-9919-07	MSI-9920-07
2	8	MSI-9917-08	MSI-9918-08	MSI-9919-08	MSI-9920-08
3	9	MSI-9917-09	MSI-9918-09	MSI-9919-09	MSI-9920-09
4	10	MSI-9917-10	MSI-9918-10	MSI-9919-10	MSI-9920-10
5	12	MSI-9917-12	MSI-9918-12	MSI-9919-12	MSI-9920-12
6	14	MSI-9917-14	MSI-9918-14	MSI-9919-14	MSI-9920-14
7	15	MSI-9917-15	MSI-9918-15	MSI-9919-15	MSI-9920-15
8	16	MSI-9917-16	MSI-9918-16	MSI-9919-16	MSI-9920-16
9	19	MSI-9917-19	MSI-9918-19	MSI-9919-19	MSI-9920-19
10	20	MSI-9917-20	MSI-9918-20	MSI-9919-20	MSI-9920-20
11	21	MSI-9917-21	MSI-9918-21	MSI-9919-21	MSI-9920-21
12	23	MSI-9917-23	MSI-9918-23	MSI-9919-23	MSI-9920-23
13	25	MSI-9917-25	MSI-9918-25	MSI-9919-25	MSI-9920-25
14	30	MSI-9917-30	MSI-9918-30	MSI-9919-30	MSI-9920-30
15	35	MSI-9917-35	MSI-9918-35	MSI-9919-35	MSI-9920-35
16	40	MSI-9917-40	MSI-9918-40	MSI-9919-40	MSI-9920-40



Fig. 00
5 mm
31 cm/12 1/4"



Fig. 0
6 mm
31 cm/12 1/4"



Fig. 1
7 mm
31 cm/12 1/4"



Fig. 2
8 mm
31 cm/12 1/4"



Fig. 3
9 mm
31 cm/12 1/4"



Fig. 4
11 mm
31 cm/12 1/4"



Fig. 5
12 mm
31 cm/12 1/4"



Fig. 6
14 mm
31 cm/12 1/4"



Fig. 7
15 mm
31 cm/12 1/4"



Fig. 8
16 mm
31 cm/12 1/4"



Fig. 9
19 mm
31 cm/12 1/4"



Fig. 10
20 mm
31 cm/12 1/4"



Fig. 11
21 mm
31 cm/12 1/4"



Fig. 12
23 mm
31 cm/12 1/4"



Fig. 13
25 mm
34 cm/13 1/2"



Fig. 14
30 mm
35 cm/12 3/4"



Fig. 15
35 mm
35 cm/12 3/4"



Fig. 16
40 mm
35 cm/12 3/4"